

# Proven™ Revision Knee System

Product Specifications



# Proven™ Revision Knee System

## Product Specifications

StelKast's patented Proven Revision Knee System features unsurpassed stem offset positioning.

### Revision Femur

- Adjustable A-P Stem Location aids in precise component fit and sizing.
- Multiple Valgus Angles – 5 & 7 Degrees Standard
- Mechanically Locked Distal & Posterior augmentation

### Revision Tibial Tray

- Adjustable Tibial Stem Offsets – 0, 2 & 4mm ensures proper fit and bone coverage
- Proximal Offset places Stem Offset Close to Tibial Plateau
- Mechanically Locked Block Augmentation
- Tri-Slot Titanium Cementless Stems
- Polished Cobalt-Chrome Cemented Stems

Proven Revision Knee System		
Catalog #	Description	Offerings
SC2236	Revision Femur	Symmetrical - Sizes 1 thru 5
SC2957		Left - Sizes 1 thru 5
SC2958		Right - Sizes 1 thru 5
SC2193	Valgus Bushing	5 & 7 Degrees Standard (0 & 3 Degrees optional)
SC2316	Distal Femoral Augment	Sizes 1 thru 5, 5 & 10mm Thickness
SC2317	Posterior Femoral Augment	Sizes 1 thru 5, 5 & 10mm Thickness
SC3029	Femoral & Tibial Stem	Cemented - 10-16mm Diameter, 40, 80 & 120mm Length
SC2196		Cementless - 10-24mm Diameter, 80 & 120mm Length
SC2241	Revision Tibial Tray	1T1F, 1T2F, 2T2F, 2T3F, 3T2F, 3T3F, 3T4F, 4T3F, 4T4F, 4T5F, 5T4F, 5T5F
SC2240	Revision Tibial Bushing	0, 2 & 4mm Offset
SC2287	½ Tibial Block (w/screws)	Sizes 1 thru 5, 5 & 10mm Thickness
SC2265	Revision Insert	Sizes 1 thru 5, 10, 12, 14, 17, 20, 23, 26 & 30mm Thickness
SC2360	Modular Femur	Sizes 1 thru 6
SC2291	Modular Tibial Tray	1T1F, 1T2F, 2T2F, 2T3F, 3T2F, 3T3F, 3T4F, 4T3F, 4T4F, 4T5F, 5T4F, 5T5F
SC2295	½ Tibial Block (w/o screws)	Sizes 1 thru 5, 5 & 10mm Thickness
SC2294	Full Tibial Wedge (w/o screws)	Sizes 1 thru 5

[www.stelkast.com](http://www.stelkast.com)



200 Hidden Valley Road | McMurray, PA 15317  
phone: (724) 941-6368 | fax: (724) 941-5987

©2008 StelKast. All rights reserved.